


DIMENSIONS FOR MATING SHAFT DETAIL PER DES-0605-M





WHEN DESIGNING SHAFT, ALLOW FOR .015 COMPRESSION BETWEEN STEP IN SHAFT AND HANDLE

REVISIONS - 601-0157

REV	DESCRIPTION	DATE/BY	DCN NO.


FOR REFERENCE ONLY



OR DENOTES INSPECTION FEATURE



THIRD ANGLE VIEW



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.

MATERIAL 6061-T6 AL, SILICONE SST END CAPS		UNSPECIFIED TOLERANCES X +/- .020 ANGLE +/- .1° XX +/- .010 MACHINE XXX +/- .005 FINISH XXXX +/- .0005		5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800	
HEAT TREAT FINISH		CONCENTRICITY .005 T.I.R. SQUARENESS .001 PER. IN.		TITLE HDL ASY, LG T F1/2-13	
DO NOT SCALE DRAWING		DRAWN HHS DATE 11/19/04		SIZE A SUPERSEDED NO.	
REMOVE MILL BURRS BREAK SHARP CORNERS		APPROVED DATE		DRAWING NO. 601-0157 REV. 00	
PART TO BE FREE OF DIRT, OIL, AND CONTAMINANTS		SCALE 1:1		ILLUSTRATION SHEET 1 OF 1	