



FOR REFERENCE ONLY


REVISIONS - BP301001-___(cc)

REV	DESCRIPTION	DATE/BY	ECN NO.



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION. THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY, WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.



BP301001-___(cc)	MATERIAL	UNSPECIFIED TOLERANCES		 beere PRECISION medical INSTRUMENTS		5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800	
	HEAT TREAT	.X ±.020 .XX ±.010 .XXX ±.005 .XXXX ±.0005	ANGLE ±1° MACHINE FINISH 32			TITLE B-AWL, PLM IMP DMD PT .146D 9.49"	
	FINISH	CONCENTRICITY .005 T.I.R. SQUARENESS .005 T.I.R.		SIZE B	SUPERSEDED NO.		
	REMOVE ALL BURRS		BREAK SHARP CORNERS		SCALE 1:1	ILLUSTRATION	
PART TO BE FREE OF DIRT, OIL AND CONTAMINANTS		APPROVED	DATE	DRAWN ASF DATE 06/20/00		DO NOT SCALE DRAWING	