


| REVISIONS | | | |
|-----------|-------------|---------|---------|
| REV | DESCRIPTION | DATE/BY | ECN NO. |
| | | | |




| STYLE | MATERIAL (BASE) | FINISH | TREATMENT (SURFACE) |
|-------|-----------------|----------------|-------------------------|
| SS | 17-4 13-8 | SATIN BRUSH | PASSIVATE PER MP-001 |
| TI | 6Al-4V | GLASS BEAD | PTFE |
| CR | 6Al-4V | GLASS BEAD | CHROME PLATE |
| NT | 6Al-4V | NONE | NONE |

NOTE:
LENGTHS CHANGE WITH POSITION OF HANDLES

SECTION A-A

THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS, INC. BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS, INC.



| | | | | |
|---|------------|--|--|--|
| MATERIAL | SEE CHART | UNSPECIFIED TOLERANCES | | BEERE PRECISION MEDICAL INSTRUMENTS, INC. 2503 LINCOLNWOOD CT., RACINE, WI 53403 (262) 554-8505 |
| | HEAT TREAT | N/A | .X +/- .020 ANGLE +/- .1° .XX +/- .010 MACHINE FINISH 32 .XXX +/- .005 .XXXX +/- .0005 | |
| FINISH | SEE CHART | CONCENTRICITY .005 T.I.R. SQUARENESS .001 PER IN. DO NOT SCALE DRAWING | TITLE PARALLEL ACTION COMPRESSOR WITH INTERCHANGEABLE HANDLES AND STANDARD BLADE SET | |
| DRAWN TJM | | DATE 8/29/01 | SIZE B | SUPERSEDED NO. |
| APPROVED | | DATE | DRAWING NUMBER PA312103- | REV. N/C |
| REMOVE ALL BURRS BREAK SHARP CORNERS PART TO BE FREE OF DIRT, OIL & CONTAMINANTS | | SCALE 1:1 | ILLUSTRATION | SHEET 1 OF 1 |