


REVISIONS			
REV	DESCRIPTION	DATE/BY	ECN NO.




NOTE:
LENGTH CHANGES WITH POSITION OF HANDLES

SECTION A-A

STYLE	MATERIAL (BASE)	FINISH	TREATMENT (SURFACE)
SS	17-4 13-8	SATIN BRUSH	PASSIVATE PER MP.001
Ti	6Al-4V	GLASS BEAD	PTFE
CR	6Al-4V	GLASS BEAD	CHROME PLATE
NT	6Al-4V	NONE	NONE

THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS, INC. BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS, INC.



MATERIAL SEE CHART	UNSPECIFIED TOLERANCES		BEERE PRECISION MEDICAL INSTRUMENTS, INC. 2503 LINCOLNWOOD CT., RACINE, WI 53403 (262) 554-8505	
	.X +/- .020	ANGLE +/- .1°	TITLE PARALLEL ACTION DISTRACTOR WITH INTERCHANGEABLE HANDLES AND EXTENDED BLADES	
.XX +/- .010	MACHINE FINISH	32		
HEAT TREAT N/A	.XXX +/- .005	CONCENTRICITY .005 T.I.R.		SIZE B
FINISH SEE CHART	.XXXX +/- .0005	SQUARENESS .001 PER IN.		SUPERSEDED NO.
	DO NOT SCALE DRAWING		DRAWING NUMBER PA312223- <u> </u>	
DRAWN TJM		DATE 8/29/01	SCALE 1:1	REV. N/C
APPROVED		DATE	ILLUSTRATION	
REMOVE ALL BURRS BREAK SHARP CORNERS		PART TO BE FREE OF DIRT, OIL & CONTAMINANTS		SHEET 1 OF 1