


FOR REFERENCE ONLY

REVISIONS - PP101112

REV	DESCRIPTION	DATE/BY	ECN NO.
C			



MATERIAL	UNSPECIFIED TOLERANCES .X +/- .020 ANGLE +/- 1° .XX +/- .010 MACHINE FINISH 32 .XXX +/- .005 .XXXX +/- .0005	BEERE PRECISION MEDICAL INSTRUMENTS 5307 95th AVENUE KENOSHA, WI 53144 (262) 657-2800
HEAT TREAT	CONCENTRICITY .005 T. I. R. SQUARENESS .001 PER IN.	TITLE PIN PULLER ASSEMBLY
FINISH	DO NOT SCALE DRAWING	SIZE SUPERSEDED NO. DWG NO. PP101112 REV C
REMOVE ALL BURRS	BREAK SHARP CORNERS	SCALE 1 : 1 ILLUSTRATION SHEET 1 OF 1
PART TO BE FREE OF DIRT, OIL, AND CONTAMINANTS	APPROVED DATE 08/19/94	

THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS, INC. BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS, INC.