


FOR REFERENCE ONLY

REVISIONS			
REV	DESCRIPTION	DATE/BY	ECO NO.
A	UPDATED ASSEMBLY WITH REV A OF P/N 601-0027	4/27/99 JMW	1346
B	UPDATED ASSEMBLY WITH REV B OF P/N 400-0001	10/5/99 ASF	1358



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS, INC. BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS, INC.

RC401000-	MATERIAL	UNSPECIFIED TOLERANCES		BEERE PRECISION MEDICAL INSTRUMENTS, INC. 2503 LINCOLNWOOD CT., RACINE WI 53403 (262) 554-8505	
	HEAT TREAT	.X +/- .020	ANGLE +/- 1°	TITLE HAND HELD ROD CUTTER, STANDARD FOR 1/4", 5.5mm, & 3/16" DIAMETER RODS	
	FINISH	.XX +/- .010	MACHINE FINISH 32		
		.XXX +/- .005	CONCENTRICITY .005 T.I.R.	SIZE B	
	.XXXX +/- .0005	SQUARENESS .001 PER IN.	DO NOT SCALE DRAWING		DATE 06/10/98
			DRAWN JAR		APPROVED
			REMOVE ALL BURRS		DATE
			BREAK SHARP CORNERS		SCALE 1:2
			PART TO BE FREE OF DIRT, OIL, AND CONTAMINANTS		ILLUSTRATION SHEET 1 OF 1