



REVISIONS - SC311xxx-___(CC)

REV	DESCRIPTION	DATE/BY	ECN NO.

FOR REFERENCE ONLY



PRODUCT SPECIFICATIONS TO BE FILLED IN AT TIME OF ORDER:


COMPLETE PART NO. _____

TORQUE SETTINGS: _____ ±10% (35 LB-IN MAX)


COLOR OF HANDLE: _____

APPROVAL/DATE: _____


OR DENOTES INSPECTION FEATURE



THIRD ANGLE VIEW



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION. THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.

MATERIAL		UNSPECIFIED TOLERANCES				5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800			
HEAT TREAT		.X +/- .020 ANGLE +/- .1° .XX +/- .010 MACHINE .XXX +/- .005 .XXXX +/- .0005 FINISH 32		TITLE T-LMT, SD, MINI 3/16" FEM SQ DR		SIZE A		SUPERSEDED NO.	
FINISH		CONCENTRICITY .005 T.I.R. SQUARENESS .001 PER. IN.		DO NOT SCALE DRAWING		DRAWING NO. SC311xxx-___(CC)		REV. E	
REMOVEMILL BURRS		BREAK SHARP CORNERS		DRAWN MTG DATE 04/15/96		SCALE 1.5:1		SHEET 1 OF 1	
PART TO BE FREE OF DIRT, OIL, AND CONTAMINANTS		APPROVED		DATE		ILLUSTRATION			