




FOR REFERENCE ONLY

REVISIONS - SP101316

REV	DESCRIPTION	DATE/BY	DCN NO.



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION. THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY, WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.

SP101316	MATERIAL	316 SST		UNSPECIFIED TOLERANCES		 beere PRECISION medical INSTRUMENTS 5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800					
	HEAT TREAT			X ±.020 XX ±.010 XXX ±.005 XXXX ±.0005	ANGLE ±1° MACHINE FINISH 32 ✓						
	FINISH			CONCENTRICITY .005 T.I.R. SQUARENESS .005 T.I.R.		TITLE SNDR, STD STR 10.75"					
	REMOVE ALL BURRS	BREAK SHARP CORNERS		DO NOT SCALE DRAWING							
PART TO BE FREE OF DIRT, OIL AND CONTAMINANTS	DRAWN	JMW	DATE	12/11/01	SIZE	B	SUPERSEDED NO.	DRAWING NO.	SP101316	REV	B
	APPROVED		DATE		SCALE	1:1	ILLUSTRATION	SHEET	1 OF 1		