





FOR REFERENCE ONLY

REVISIONS - SP102316

REV	DESCRIPTION	DATE/BY	DCN NO.



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION. THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY, WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.

SP102316	MATERIAL	316 SST		UNSPECIFIED TOLERANCES		 5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800	
	HEAT TREAT			X ±.020 XX ±.010 XXX ±.005 XXXX ±.0005	ANGLE ±1° MACHINE FINISH 32 ✓		
	FINISH			CONCENTRICITY .005 T.I.R. SQUARENESS .005 T.I.R.		TITLE SNDR, STD CRV	
	REMOVE ALL BURRS		BREAK SHARP CORNERS		DO NOT SCALE DRAWING		
PART TO BE FREE OF DIRT, OIL AND CONTAMINANTS		DRAWN	DAM	DATE	05/10/05	SIZE	B
		APPROVED		DATE		SUPERSEDED NO.	
						DRAWING NO.	SP102316
						SCALE	1:1
						ILLUSTRATION	
						SHEET	1 OF 1
						REV	B