



# FOR REFERENCE ONLY

REVISIONS - ST301\_\_\_\_(length mm)-\_\_(cc)



REV	DESCRIPTION	DATE/BY	ECN NO.



- LASER MARK DATA AND SEQUENCE (left to right) 180° FROM LENGTH MARKINGS  
 CHECK DATA REQUIRED (\*REQUIRED FOR ALL PARTS)
- \_\_\_\_ PART NUMBER \_\_\_\_\_
  - \_\_\_\_ REVISION \_\_\_\_\_
  - \_\_\_\_ CUSTOMER LOGO? YES / NO (circle one)
  - \_\_\_\_ CE MARK? YES / NO (circle one) REG. NUMBER \_\_\_\_\_
  - \_\_\_\_ LOT NUMBER \_\_\_\_\_
  - \_\_\_\_ \* "FOR SINGLE USE ONLY - DO NOT IMPLANT OR CUT - PATENTED"
  - \_\_\_\_ OTHER \_\_\_\_\_

APPROVED BY (INITIAL & DATE): \_\_\_\_\_

OR DENOTES INSPECTION FEATURE



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION. THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.

MATERIAL SEE BOM		UNSPECIFIED TOLERANCES .X +/- .020 ANGLE +/- 1° .XX +/- .010 MACHINE 32 .XXX +/- .005 .XXXX +/- .0005 FINISH		 5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800	
HEAT TREAT		CONCENTRICITY .005 T.I.R. SQUARENESS .001 PER. IN.		TITLE T-ROD, 3.0mm SIL COVERED ANDZ ALUM	
FINISH		DO NOT SCALE DRAWING		DRAWING NO. ST301____(length mm)-__(cc)	
REMOVEMILL BURRS	BREAK SHARP CORNERS	DRAWN DRP	DATE 03/01/02	SIZE A	REV B
PART TO BE FREE OF DIRT, OIL, AND CONTAMINANTS		APPROVED	DATE	SCALE 1:1	SHEET 1 OF 1