


FOR REFERENCE ONLY

REVISIONS - ST601____(length mm)-____(cc)

REV	DESCRIPTION	DATE/BY	ECN NO.





PART NUMBER SHOWN FOR REFERENCE ONLY, FOR ACTUAL PART NUMBER SEE INSTRUCTIONS BELOW

LASER MARK DATA AND SEQUENCE (left to right) 180° FROM LENGTH MARKINGS CHECK DATA REQUIRED (*REQUIRED FOR ALL PARTS)

- ___ PART NUMBER _____
- ___ REVISION _____
- ___ CUSTOMER LOGO? YES / NO (circle one)
- ___ CE MARK? YES / NO (circle one) REG. NUMBER _____
- ___ LOT NUMBER _____
- ___ * "FOR SINGLE USE ONLY - DO NOT IMPLANT OR CUT - PATENTED"
- ___ OTHER _____

APPROVED BY (INITIAL & DATE): _____

OR DENOTES INSPECTION FEATURE



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS.

MATERIAL SEE B.O.M		UNSPECIFIED TOLERANCES X +/- .020 ANGLE +/- .1° XX +/- .010 MACHINE XXX +/- .005 FINISH XXXX +/- .0005				5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800	
HEAT TREAT		CONCENTRICITY .005 T.I.R. SQUARENESS .001 PER. IN.		TITLE T-ROD, 6.0mm SIL COVERED ANDZ ALUM		DRAWING NO. ST601____(length mm)-____(cc)	
FINISH		DO NOT SCALE DRAWING		SIZE A		REV C	
REMOVEMILL BURRS	BREAK SHARP CORNERS	DRAWN JAR	DATE 08/20/99	SUPERSEDED NO.		SHEET 1 OF 1	
PART TO BE FREE OF DIRT, OIL, AND CONTAMINANTS		APPROVED	DATE	SCALE 1:1		ILLUSTRATION	