






FOR REFERENCE ONLY

REVISIONS - TC104000


REV	DESCRIPTION	DATE/BY	ECN NO.



LASERMARK "PART NUMBER", "SHOP ORDER NUMBER", AND "NAME" IN .10 HIGH CHARACTERS. ORIENTATE WITH FLATS AND BALL AS SHOWN. "USE BOLD TEXT".



THIS DOCUMENT INCLUDING ALL DATA, SPECIFICATIONS, AND INFORMATION CONTAINED WITHIN IS CONFIDENTIAL AND PROPRIETARY TO BEERE PRECISION MEDICAL INSTRUMENTS, INC. BY RETENTION, THE RECIPIENT AGREES NOT TO REPRODUCE, DISCLOSE, OR MAKE USE OF ITS CONTENTS IN ANY WAY WITHOUT PRIOR WRITTEN PERMISSION OF AN OFFICER OF BEERE PRECISION MEDICAL INSTRUMENTS, INC.

TC104000	MATERIAL	UNSPECIFIED TOLERANCES		 beere PRECISION medical INSTRUMENTS 5307 95TH AVE KENOSHA, WI 53144 (262) 657-2800			
	HEAT TREAT	.X ±.020 .XX ±.010 .XXX ±.005 .XXXX ±.0005	ANGLE ±1° MACHINE FINISH 32/				
	FINISH	CONCENTRICITY .005 T.I.R. SQUARENESS .005 T.I.R.		TITLE DVR ASY, C-TORQ			
	REMOVE ALL BURRS BREAK SHARP CORNERS PART TO BE FREE OF DIRT, OIL AND CONTAMINANTS	DO NOT SCALE DRAWING	DRAWN JAR	DATE 12/27/01	SIZE B	SUPERSEDED NO.	
				SCALE 1:1	DRAWING NO. TC104000	REV N/C	
				SHEET 1 OF 1		ILLUSTRATION	